


Fab and Fab-PEG(WCX)



Column: BioCore WCX, 10 μ m
Dimension: 4.6 \times 150 mm
Mobile Phase: A) 20 mM MES, pH5.5
B) 300 mM NaCl in 20 mM MES, pH5.5
Gradient:

| t(min) | %A | %B |
|--------|-----|-----|
| -10 | 100 | 0 |
| 0 | 100 | 0 |
| 20 | 60 | 40 |
| 20.1 | 0 | 100 |
| 25 | 0 | 100 |

Flow Rate: 1.0 mL/min
Temperature: 30 $^{\circ}$ C
Injection: 2 μ L
Detection: UV 280 nm
Sample: Fab-PEG (3 mg/mL)
Fab (5 mg/mL)