


ADC与Pro蛋白



Column: BioCore HIC-Butyl, 5 μ m
Dimension: 4.6 \times 100 mm
Mobile Phase: A) 2.0 M $(\text{NH}_4)_2\text{SO}_4$ in 100 mM phosphate buffer, pH7.0
B) 100 mM phosphate buffer, pH7.0
C) Isopropanol

Gradient:	t(min)	%A	%B	%C
	0	50	50	0
	1	50	50	0
	15	0	75	25
	20	0	75	0
	20.1	50	50	0
	35	50	50	0

Flow Rate: 1.0 mL/min
Temperature: 30 $^\circ\text{C}$
Injection: 20 μL
Detection: UV 220 nm
Samples: ADC (0.46 mg/mL)
Pro蛋白 (1.8 mg/mL)