


双抗



Columns:

Blue: BioCore HIC-Phenyl, 5 μ m

Black: BioCore HIC-Butyl, 5 μ m

Dimension: 4.6 \times 100 mm

Mobile Phase: A) 7.5/92.5 v/v Isopropanol/2.0 M $(\text{NH}_4)_2\text{SO}_4$ in 100 mM phosphate buffer, pH7.0

B) 7.5/92.5 v/v Isopropanol/100 mM phosphate buffer, pH7.0

Gradient:

t (min)	%A	%B
0	100	0
1	100	0
20	0	100
20.1	100	0
26	100	0

Flow Rate: 1.0 mL/min

Temperature: 30 $^{\circ}$ C

Injection: 2 μ L

Detection: UV 280 nm

Instrument: Agilent 1290 Infinity II (Bio)

Sample: 双抗 (1.72 mg/mL)

